


Super Medium Disaster Pouch, Envelope Zipper, HDPE Tarp, 6 Handles | Size Options

Mopec is proud to partner with Salam International to bring you the best body bags available in the market.



SPECIFICATIONS

Construction

- Ends are high frequency welded
- Heavy-duty #10 zipper
- 2 lockable pulls to secure pouch
- Complies with OSHA regulation for containment of body fluids, blood, body fats, chemicals, and fluids
- Envelope Zipper
- Color: Dark blue
- 6 Handles

Materials

- HDPE specially woven coated tarpaulin
- 2" wide webbing for easy grabbing
- Chlorine free

Dimensions

- 94" L x 36" W (BE411-324-EZ)
- 94" L x 46" W (BE411-354-EZ)

Performance

- 2 way zipper so pouch can open and close at any point
- No sewing or welding on sides to increase strength and eliminate leaks
- Handles secured by heavy-duty rivets
- Handle straps run under the pouch for additional strength and support
- 100% leakproof
- Recommended for cremation

Options & Accessories

- Sold by the case of 15 (BE411-324-EZ) or 5 (BE411-354-EZ)
- Available in adult (BE411-324-EZ) and bariatric (BE411-354-EZ) sizes.

Options

Part No.	Dimensions	Zipper Style	Color	Handles
BE411-324-EZ	94" L x 36" W	Envelope	Dark Blue	6 Handles
BE411-354-EZ	94" L x 46" W	Envelope	Dark Blue	6 Handles

SALAM
SALAM INTERNATIONAL, INC.

AVAILABLE FROM
Mopec