

Pedestal Autopsy Table, L-Shaped, Right or Left Facing, with Optional Elevation

1036-19AR



1036-19L



SPECIFICATIONS

Construction

- Heliarc-welded seams and joints, ground and polished to match adjoining surfaces
- Large single sink basin
- Mixing faucet with gooseneck spout and wrist blade handles
- (4) Fully perforated grid plates
- Heavy duty waste disposal (1/2 HP)
- Deluxe hydro-aspirator with reverse flow, cold water control valve with 10ft of hose
- Spray hose assembly with cold water control valve, nozzle with 10ft of hose

Materials

- 304 stainless steel construction with a polished No. 4 finish

Dimensions

- Large side: 104"L x 32"W x 35"H
- Small side: 86 1/16"L x 24"W x 35"H
- Elevation Option: 32"-44"H
- Sink: 14"L x 12"W x 8"D

Utilities / Power

- (1) 115v/60hz/1ph/15a circuit for station and (1) 115v/60hz/1ph/15a circuit for disposer
- Exhaust recommendation: 450-500 CFM

Capacity

- 750 lbs

Performance

- GFCI protected receptacles with waterproof covers
- Back-flow protection for all plumbing fixtures to prevent back-flow of contaminated water
- Rapid and positive drainage
- (2) Access panels on pedestal

Options & Accessories

- Elevating Units: Select 1036-19AL or 1036-19AR
- For a list of available accessories, contact your representative

Related/Alternate/Obsolete Product Model Numbers

Model No.	Description	Status	Notes
CE600	Pedestal Autopsy Table, L-Shaped, Right	Replaced	Use 1036-19R
CE650	Pedestal Autopsy Table, L-Shaped, Elevating, Right	Replaced	Use 1036-19AR
CE700	Pedestal Autopsy Table, L-Shaped, Left	Replaced	Use 1036-19L
CE750	Pedestal Autopsy Table, L-Shaped, Elevating, Left	Replaced	Use 1036-19AL
1036-19-IHA	Pedestal Autopsy Table, L-Shaped, Elevating	Obsolete	
1036-19-SL	Pedestal Autopsy Table, L-Shaped, Side Loading	Obsolete	

Variations

Model No.	Description	Elevation
1036-19L	Pedestal Autopsy Table, L-Shaped, Left	Non-Elevating
1036-19R	Pedestal Autopsy Table, L-Shaped, Right	Non-Elevating
1036-19AL	Pedestal Autopsy Table, L-Shaped, Elevating, Left	Elevating
1036-19AR	Pedestal Autopsy Table, L-Shaped, Elevating, Right	Elevating

MOPEC
GROUP

mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI