

Dissection Table with Hand-Crank Elevation - 600023HC

The 600023HC Dissection Table is designed to optimize comfort, safety, and efficiency in anatomy labs. Its adjustable height ensures ergonomic support for users, while the versatile and mobile design, enhances workflow and reflects a commitment to innovative lab solutions.



SPECIFICATIONS

Construction

- Hand-crank elevation
- Recessed tabletop, creased towards drain for added fluid flow
- 6" casters, all with brake mechanisms
- Accessory holder (support sockets) are located on all (4) corners to accommodate optional accessories
- Full-length utility shelf

Materials

- 304 stainless steel construction with a polished No. 4 finish

Dimensions

- 80" L x 30" W x 32"-44" H

Capacity

- 350 Lbs

Options & Accessories

- For a list of accessories, speak with your Mopec Group representative

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI