


Short Shelf for Dissection Table



SPECIFICATIONS

Construction

- Spans 3/4 length of the underside of your dissection table
- Meant for use with AX107: bucket holder

Materials

- 304 stainless steel construction with a polished No. 4 finish

Dimensions

- 66.3" L x 23.3" W x 1" H
- Usable area: 66.2" L x 19.6" W

Compatibility

- HA100
- HA210
- HA400
- HB1400
- HB1410
- HB1500
- HB1510

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI