


# Ensemble Pegboard Marking Dye Holder



## SPECIFICATIONS

### Construction

- Designed for the Mopec Maestro's Ensemble Pegboard
- Fits 9 large dye bottles less than 2.0" in diameter, and marking dye applicator sticks

### Materials

- 304 stainless steel construction with a polished No. 4 finish

### Dimensions

- 12.3" W x 3.0" H x 5.3" D
- 9 dye holes: 2.0" diameter
- 2 stick holes: .75" diameter

### Compatibility

- Fits all Mopec Pegboards, including AX031R and AX031L pegboard wings
- Can also be attached to the front of grossing stations, above the vent grille

**MOPEC**  
GROUP

**Mopec**  
ELEVATING PATHOLOGY

**MORTECH**  
MANUFACTURING

[www.mopecgroup.com](http://www.mopecgroup.com)

800.362.8491

800 Tech Row, Madison Heights, MI