

Ensemble Pegboard Bins



AX008



AX010

SPECIFICATIONS

Construction

- Designed for the Mopec Maestro's Ensemble Pegboard

Materials

- 304 stainless steel construction with a polished No. 4 finish

Dimensions

- AX008: 4.0" L x 4.2" W x 4.0" H
- AX010: 5.0" L x 5.2" W x 5.0" H

Capacity

- AX008: 10 lbs
- AX010: 12 lbs

Compatibility

- Fits all Mopec Pegboards, including AX031R and AX031L pegboard wings
- Can also be attached to the front of grossing stations, above the vent grille

Variations

Model No.	Description	Dimensions	Capacity
AX008	Pegboard Bin, Small	4" L x 4.2" W x 4" H	10 lbs
AX010	Pegboard Bin, Large	5" L x 5.2" W x 5" H	12 lbs

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI