


# Ensemble Pegboard C-Fold Paper Towel & Glove Bins




AX011



AX013



AX012



AX014

## SPECIFICATIONS

### Construction

- Designed for the Mopec Maestro's Ensemble Pegboard
- Designed to hold boxes of gloves or stacks of C-fold paper towels.

### Materials

- 304 stainless steel construction with a polished No. 4 finish

### Dimensions

- AX011: 10.8" L x 4.2" W x 5.0" H
- AX012: 10.8" L x 4.2" W x 10.0" H
- AX013: 15.8" L x 4.0" W x 10.8" H
- AX014: 20.8" L x 4.0" W x 10.8" H

### Capacity

- AX011: 10 lbs
- AX012: 15 lbs
- AX013: 16 lbs
- AX014: 18 lbs

### Compatibility

- Fits all Mopec Pegboards, including AX031R and AX031L pegboard wings

## Variations

| Model No. | Description                         | Dimensions                 |
|-----------|-------------------------------------|----------------------------|
| AX011     | Pegboard C-Fold / Glove Bin, Single | 10.8" L x 4.2" W x 5" H    |
| AX012     | Pegboard C-Fold / Glove Bin, Double | 10.8" L x 4.2" W x 10" H   |
| AX013     | Pegboard C-Fold / Glove Bin, Triple | 15.8" L x 4.0" W x 10.8" H |
| AX014     | Pegboard C-Fold / Glove Bin, Quad   | 20.8" L x 4.0" W x 10.8" H |

**MOPEC**  
GROUP

**mopec**  
ELEVATING PATHOLOGY

**MORTECH**  
MANUFACTURING

[www.mopecgroup.com](http://www.mopecgroup.com)

800.362.8491

800 Tech Row, Madison Heights, MI