


Mopec Maestro Encore PathCamHR Factory Kit



SPECIFICATIONS

Construction

- (1) 15a receptacle, 125 vac, 35mm din rail mount to supply camera power
- (1) Panel-Mount USB-A 2.0 Coupler with cover plate
- Includes all necessary hardware to complete installation

Materials

- Bracket: 304 stainless steel construction

Dimensions

- Mounting plate: 5.8" L x 8.8" W x 4.6" H
- Receptacle: 5.7" L x 2.1" W x 2.3" H
- Panel-Mount USB with cover plate: 1.7" L x 2.0" W x 4.3" H

Compatibility

- Included on all Maestro Encore platforms
- AX265HD camera
- AX266HR camera

MOPEC
GROUP

mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI