

Post-Style Leg Holder for Dissection Table



SPECIFICATIONS

Construction

- Includes 2 leg holders and support posts
- Built to support each leg while working
- Leg holders can be used on any of the 4 corners of the cart

Materials

- 304 stainless steel construction with a polished No. 4 finish

Dimensions

- 4.3" L x 5.6" W x 24.9" H

Compatibility

- HA100
- HA200
- HA210
- HA400
- HB1400
- HB1410
- HB1500
- HB1510
- HB1550

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI