

# Ensemble Pegboard Magnetic Toolbar



## SPECIFICATIONS

### Construction

- Designed for the Mopec Maestro's Ensemble Pegboard
- Powerful magnet for attaching all your tool and instruments

### Capacity

- 5 lbs

### Materials

- 304 stainless steel construction with a polished No. 4 finish
- Rare-earth magnetic surface

### Dimensions

- Overall: 14.8" L x .8" W x 1.3" H
- Magnet: 13.3" L x 1" H

### Compatibility

- Fits all Mopec Pegboards, including AX031R and AX031L pegboard wings
- Can also be attached to the front of grossing stations, above the vent grille

**MOPEC**  
GROUP

**Mopec**  
ELEVATING PATHOLOGY

**MORTECH**  
MANUFACTURING

[www.mopecgroup.com](http://www.mopecgroup.com)

800.362.8491

800 Tech Row, Madison Heights, MI