





CUTTING BOARD IS TO BE FREE FROM SCRATCHES, STAINS, DISCOLORATION, & DEFORMATION

REV	DESCRIPTION	DATE	ECN #	APPR
A	INITIAL RELEASE	00-XXX-19	XXXXX	MM




THIS DOCUMENT CONTAINS PROPRIETARY INFORMATION OF MOPEC, AND MAY NOT BE COPIED IN WHOLE OR PART, AND NO USE SHALL BE MADE OF IT OTHER THAN FOR WHICH IT HAS BEEN SUBMITTED WITHOUT THE WRITTEN PERMISSION OF A MOPEC MANAGER MOPEC 2021	UNLESS OTHERWISE SPECIFIED ALL DIMENSIONS IN INCHES TOLERANCES 0.0 ± 0.1 in 0.00 ± 0.05 in 0.000 ± 0.010 in 0.0000 ± 0.0010 in ANGLES ±0.5° BREAK ALL CORNERS 0.063 X 45° ALL FILLET RADII R0.063 /0.125in INSPECT OVAL DIMENSIONS	800 Tech Row, Madison Heights, MI 48071	DRAWN BY A. SALAME DATE 4/11/2022	PART NUMBER AX064-048-A
			CHECKED BY DATE 4/11/2022	DESCRIPTION MAESTRO, DISSECTION BOARD, WHITE, 48"
DRAWN TO 3RD ANGLE PROJECTION			APPROVED BY DATE: M.McCLAIN 4/11/2022	SIZE A2 FILE NAME AX064-048-A.dft REV A
SCALE 1:2			DO NOT SCALE IF IN DOUBT ASK	SHEET 1 OF 2

AX064-048-A.dft

CUTTING BOARD IS TO BE FREE FROM SCRATCHES, STAINS, DISCOLORATION, & DEFORMATION

REV	DESCRIPTION	DATE	ECN #	APPR
A	INITIAL RELEASE	00-XXX-19	XXXXX	MM



Item Number	Quantity	File Name	Title	Material	Revision
1	4	PB0609.par	TAPPING INSERT, 1/4"-20 INTERNAL, BRASS	McMASTER-CARR # 90016A029	A
2	1	S7541-A.par	DISSECTION BOARD, 15" x 18.5", WHITE	3/4" HDPE - NATURAL, TEXTURED	A

THIS DOCUMENT CONTAINS PROPRIETARY INFORMATION OF MOPEC, AND NO USE SHALL BE MADE OF IT OTHER THAN FOR WHICH IT HAS BEEN SUBMITTED WITHOUT THE WRITTEN PERMISSION OF A MOPEC MANAGER
MOPEC 2021

UNLESS OTHERWISE SPECIFIED ALL DIMENSIONS IN INCHES
TOLERANCES 0.0 ± 0.1 in
0.00 ± 0.05 in
0.000 ± 0.010 in
0.0000 ± 0.0010 in
ANGLES ±0.5°
BREAK ALL CORNERS 0.063 X 45°
ALL FILLET RADII R0.063 /0.125in
INSPECT OVAL DIMENSIONS

Mopec
800 Tech Row,
Madison Heights, MI
48071
DRAWN TO 3RD ANGLE PROJECTION

DRAWN BY: A. SALAME
DATE: 4/11/2022
CHECKED BY: ---
DATE: 4/11/2022
APPROVED BY: M. McCLAIN
DATE: 4/11/2022

PART NUMBER: AX064-048-A
DESCRIPTION: MAESTRO, DISSECTION BOARD, WHITE, 48"
SIZE: A2 FILE NAME: AX064-048-A.dft REV: A
SCALE: 1:2.67 DO NOT SCALE IF IN DOUBT ASK SHEET 2 OF 2

AX064-048-A.dft