

Ensemble Dissection Board, Light Blue, with 3” Legs



SPECIFICATIONS

Construction

- Made to the depth of the grossing station working area and is the width of one half of the grossing station work surface

Dimensions

- AX067-060: 27” L 18.5” D X 0.8” H
- For additional sizes, contact your Mopec representative

Capacity

- 20 lbs

Materials

- HDPE, High-Density polyethylene, light blue in color
- HDPE legs with rubber feet

Compatibility

- Fits all Mopec Maestro and Maestro Encore grossing stations

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI