


REV	DESCRIPTION	DATE	AP
A	INITIAL RELEASE	26-DEC-18	M



Item Number	Quantity	File Name	Title	Material	Rev	<div>THIS DOCUMENT CONTAINS PROPRIETARY INFORMATION OF MOPEC, AND MAY NOT BE COPIED IN WHOLE OR PART, AND NO USE SHALL BE MADE OF IT OTHER THAN FOR WHICH IT HAS BEEN SUBMITTED WITHOUT THE WRITTEN PERMISSION OF A MOPEC MANAGER</div> <div>MOPEC 2024</div>	<div>UNLESS OTHERWISE SPECIFIED ALL DIMENSIONS IN INCHES TOLERANCES</div> <div>0.00 ± 0.2 in 0.00 ± 0.10 in 0.000 ± 0.005 in 0.0000 ± 0.0005 in ANGLES ±0.5°</div> <div>BREAK ALL CORNERS 0.063 X 45° ALL FILLET RADI R0.063 /0.125in INSPECT OVAL DIMENSIONS</div>	<div><div>MOPEC GROUP</div><div>MopecMORTECH</div></div> <div><div>DRAWN BY DATE: B. LAY 26-DEC-18</div><div>CHECKED BY DATE: --- 26-DEC-18</div><div>APPROVED BY DATE: M.McCLAIN 26-DEC-18</div></div> <div><div>DRAWN TO 3RD ANGLE PROJECTION</div><div></div></div>	<div>PART NUMBER AX070</div> <div>DESCRIPTION ADJ. MONITOR AND KEYBOARD ARMS</div> <div>SIZE A2 FILE NAME AX070.dft</div> <div>SCALE 1:1</div>			<div>REV A</div> <div>DO NOT SCALE IF IN DOUBT ASK</div> <div>SHEET 2 OF 2</div>	
1	1	pm0354.asm	LX KEYBOARD ARM WALL MOUNT										
2	1	pm0355-h.asm	LX WALL MOUNT LCD ARM		H								
3	2	S6836-B.psr	AX ACCESSORY MOUNT	RM0258, FLAT, AL, 1/2" X 2-1/2", 6061	B								

AX070.dft