



# Ensemble Cassette Labeling Shelf



AX073L



AX073R

## SPECIFICATIONS

### Construction

- Designed to accommodate the most popular cassette labelers

### Materials

- 304 stainless steel construction with a polished No. 4 finish

### Dimensions

- 13.1" L x 32" D x 11" H
- Usable area: 12.9" L x 30.8" D

### Performance

- Either option can be used for left-handed or right-handed Maestro and Maestro Encore grossing stations

### Capacity

- 40 lbs

### Compatibility

- Fits all Mopec Maestro and Maestro Encore grossing stations

## Variations

| Model No. | Description                    | Dimensions              |
|-----------|--------------------------------|-------------------------|
| AX073L    | Cassette Labeling Shelf, Left  | 13.1" L x 32" D x 11" H |
| AX073R    | Cassette Labeling Shelf, Right | 13.1" L x 32" D x 11" H |

**MOPEC**  
GROUP

**mopec**  
ELEVATING PATHOLOGY

**MORTECH**  
MANUFACTURING

[www.mopecgroup.com](http://www.mopecgroup.com)

800.362.8491

800 Tech Row, Madison Heights, MI