

# Ensemble Shelf-Mounted Single Gooseneck LED Task Light



## SPECIFICATIONS

### Construction

- Mounts to the shelf area of the grossing station directly in the middle of the working area

### Materials

- ABS plastic, rubber coated flexible metal arm, stainless steel hardware

### Dimensions

- 6.5" L x 2.4" W x 17.3" H
- Flexible arm: 15.3"
- Flexible arm with light head: 16.5"

### Compatibility

- Fits all Mopec Maestro and Maestro Encore grossing stations

### Options & Accessories

- For a pair of lights, please order AX082 Dual Gooseneck LED Task Lights

## NOTE:

- We recommend ordering at the time of purchase to be properly installed during manufacturing of the grossing station.

**MOPEC**  
GROUP

**mopec**  
ELEVATING PATHOLOGY

**MORTECH**  
MANUFACTURING

[www.mopecgroup.com](http://www.mopecgroup.com)

800.362.8491

800 Tech Row, Madison Heights, MI