

PathCam Hands-Free Voice Command



SPECIFICATIONS

Construction

- Designed to allow the user to program voice commands for functions of their choosing
- Includes wireless headset microphone and USB dongle connection

Performance

- Preprogrammed and tested
- Users may customize voice commands or use default commands, including “zoom in”, “zoom out”, and “snapshot”

Compatibility

- Maestro and Maestro Encore grossing stations equipped with PathCam Imaging Systems

MOPEC
GROUP

iMopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI