


Maestro Recirculating to In-House Facility Vented Conversion Kit



SPECIFICATIONS

Construction

- Designed to convert recirculating grossing stations into in-house ventilation
- Bolt-on panel with 8" duct stubs which allows the grossing station to be connected to the facility ventilation system
- Includes panel, duct tubing, and clamps

Dimensions

- 10.2" D X 4.5" H
- AX110-048: 41" L
- AX110-060: 53" L
- AX110-072: 65" L
- AX110-096: 89.1" L

Materials

- 304 stainless steel sheetmetal and hardware

Compatibility

- Fits the Mopec Maestro and Maestro Encore grossing stations of the accompanying sized unit

Variations

Model No.	Description	Dimensions
AX110-048	Recirculating to In-House Kit, 48" Units	41" L x 10.2" D x 4.5" H
AX110-060	Recirculating to In-House Kit, 60" Units	53" L x 10.2" D x 4.5" H
AX110-072	Recirculating to In-House Kit, 72" Units	65" L x 10.2" D x 4.5" H
AX110-096	Recirculating to In-House Kit, 96" Units	89.1" L x 10.2" D x 4.5" H

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI