

# Ensemble Dual Monitor Mount Arm



\*Monitors not included

## SPECIFICATIONS

### Construction

- Designed to mount two monitors into the existing track on the column opposite the sink
- Once installed, track allows arm to be moved up and down the track by loosening two large, knurled knobs, without needing any tools

### Materials

- Keyboard arm components are powder-coated aluminum and steel for the brackets

### Dimensions

- 27.5" L x 13.9" D x 17.2" H

### Performance

- Has several pivot points to articulate up and down, and side to side.
- Can be pulled out when in use, then pushed back when not in use
- Arm has 24" of extension, 180° of monitor pan, and 85° of tilt

### Capacity

- 26 lbs total, 13 lbs per monitor

### Compatibility

- Fits all Mopec Maestro and Maestro Encore grossing stations
- Safely accommodates 22" and 24" monitors (maximum width of 21.25")
- **Note:** monitors not included

**MOPEC**  
GROUP

**mopec**  
ELEVATING PATHOLOGY

**MORTECH**  
MANUFACTURING

[www.mopecgroup.com](http://www.mopecgroup.com)

800.362.8491

800 Tech Row, Madison Heights, MI