

Ensemble Magnifier with LED Light



SPECIFICATIONS

Construction

- Designed to mount on the shelf of a Maestro or Maestro Encore grossing station to help better see specimens at a 1.75x magnification

Materials

- Powder coated die cast aluminum
- Glass 1.75X diopter lens

Dimensions

- 6.8" L x 21.1" D x 14.9" H

Performance

- Magnifies 1.75x
- Light has color of 4000K, 700 lumens, is dimmable, and has a 50,000-hour lamp life
- Arm has 30" of extension and 180° of sweep

Compatibility

- Fits all Mopec Maestro and Maestro Encore grossing stations

NOTE:

- This item **must** be ordered at the time of purchase to be properly installed during manufacturing of the grossing station.

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI