

Stainless Steel Supply Cart & Podium

Designed for seamless mobility, this dual-purpose unit provides a durable, easy-to-clean workspace while keeping essential supplies within reach, making it ideal for laboratory environments. Built with Mopec's commitment to quality and innovation, this cart optimizes workflow and storage while maintaining a professional and ergonomic design.



SPECIFICATIONS

Construction

- Heliarc-welded seams and joints, ground and polished to match adjoining surfaces
- (4) 3" locking casters
- Built-in shelving

Materials

- 304 stainless steel construction with a polished No. 4 finish

Performance

- Angled writing surface with book stop and pencil ledge

Dimensions

- 24.3" L x 24" W x 48.6" H
- Top Shelf: 21.9" W x 11" H
- Middle Shelf: 21.9" W x 11" H
- Bottom Shelf: 21.9" W x 16" H

Capacity

- 80 lbs

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI