

Stainless Steel Book Holding Stations

Optimize organization and accessibility with Mopec’s Book Holding Stations, built for laboratories, pathology suites, and educational environments. Both the mobile and stationary options offer a versatile reference podium that feature a unique dual-sided book rest—perfect for teaching applications and post-mortem documentation. Its extended bottom lips securely hold thick medical publications, paperwork, or e-readers, while two storage trays provide space for essential tools and devices.



SPECIFICATIONS

Construction

- Heliarc-welded seams and joints, ground and polished to match adjoining surfaces
- (2) Angled book surfaces with book edge
- (2) Flat storage areas with 2” H walls

Materials

- 304 stainless steel construction with a polished No. 4 finish

Performance

- (2) threaded knobs for adjusting height
- BN100: (2) 3.5” swivel front casters and (2) 3.5” locking swivel rear casters
- BN200: (4) adjustable leveling feet

Dimensions

- BN100: 24” L x 20” W x 51 - 55” H
- BN200: 24” L x 20” W x 47 - 50” H
- Angled surfaces: 20” W x 12.4” H
- Flat areas: 17.6” L x 13.9” W x 2” D

Capacity

- 80 lbs

Options & Accessories

- BN100: Mobile
- BN200: Stationary

Variations

Model No.	Description	Dimensions
BN100	Book Holding Station, Mobile	24” L x 20” W x 51-55” H
BN200	Book Holding Station, Stationary	24” L x 20” W x 47-50” H

MOPEC
GROUP

mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI