

Maestro Cabinet Doors



MB1023-048



MB1023-060



MB1023-072



MB1023-096

SPECIFICATIONS

Construction

- Designed to close off the “knee space” of your Maestro or Maestro Encore grossing station
- Great for keeping formalin dispensing cubitainer, formalin collection container, or other storage out of sight
- Can be added to your station at any time

Materials

- 304 stainless steel construction with a No. 4 finish

Dimensions

- 2.3” D x 25.7” H
- MB1023-048: 22.8” L
- MB1023-060: 34.8” L
- MB1023-072: 46.8” L
- MB1023-096: 29.5” L

Quantity

- MB1023-048: 1 cabinet door
- MB1023-060: 1 pair of doors
- MB1023-072: 1 pair of doors
- MB1023-096: 2 pairs of doors

Compatibility

- Compatible with all Maestro and Maestro Encore grossing stations
- MB1023-048: 48” stations
- MB1023-060: 60” stations
- MB1023-072: 72” stations
- MB1023-096: 96” stations

Variations

Model No.	Description	Quantity	Dimensions
MB1023-048	Cabinet Door, 48” Units	1 Door	22.8” L x 2.3” D x 25.7” H
MB1023-060	Cabinet Doors, 60” Units	2 Doors	34.8” L x 2.3” D x 25.7” H
MB1023-072	Cabinet Doors, 72” Units	2 Doors	46.8” L x 2.3” D x 25.7” H
MB1023-096	Cabinet Doors, 96” Units	4 Doors	29.5” L x 2.3” D x 25.7” H

MOPEC
GROUP

mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI